




FAF1 Antibody

Product Code	CSB-PA892359LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNN5
Immunogen	Recombinant Human FAS-associated factor 1 protein (490-629AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAS-associated factor 1 (hFAF1) (UBX domain-containing protein 12) (UBX domain-containing protein 3A), FAF1, UBXD12, UBXN3A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAF1

Image




Western Blot

Positive WB detected in: HeLa whole cell lysate, HEK293 whole cell lysate, THP-1 whole cell lysate, Mouse heart tissue, Mouse brain tissue
All lanes: FAF1 antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 74, 57 kDa
Observed band size: 74 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA892359LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.