




# FAF1 Antibody

<b>Product Code</b>	CSB-PA892359LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UNN5
<b>Immunogen</b>	Recombinant Human FAS-associated factor 1 protein (490-629AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	FAS-associated factor 1 (hFAF1) (UBX domain-containing protein 12) (UBX domain-containing protein 3A), FAF1, UBXD12, UBXN3A
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	FAF1

## Image




### Western Blot

Positive WB detected in: Hela whole cell lysate, HEK293 whole cell lysate, THP-1 whole cell lysate, Mouse heart tissue, Mouse brain tissue  
All lanes: FAF1 antibody at 2.8μg/ml

### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution  
Predicted band size: 74, 57 kDa  
Observed band size: 74 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA892359LA01HU at dilution of 1:100