



FAF1 Antibody, HRP conjugated

Product Code	CSB-PA892359LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNN5
Immunogen	Recombinant Human FAS-associated factor 1 protein (490-629AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAS-associated factor 1 (hFAF1) (UBX domain-containing protein 12) (UBX domain-containing protein 3A), FAF1, UBXD12, UBXN3A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.