



BTRC Antibody, Biotin conjugated

Product Code	CSB-PA897464LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y297
Immunogen	Recombinant Human F-box/WD repeat-containing protein 1A protein (397-448AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 1A (E3RSIkappaB) (Epididymis tissue protein Li 2a) (F-box and WD repeats protein beta-TrCP) (pIkappaBalph-E3 receptor subunit), BTRC, BTRCP FBW1A FBXW1A
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	BTRC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.