



Recombinant Human Intestinal-type alkaline phosphatase (ALPI), partial

Product Code	CSB-YP001627HU1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09923
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ENPAFWNRQAAEALDAAKKLQPIQKVAKNLILFLGDGLGVPTVTATRILKGQKN GKLGPEPLAMDRFPYLALSPTYNVDRQVPDSAATATAYLCGVKANFQTIGLS AAARFNQCNTTRGNEVISVMNRAKQAGKSVGVTTRVQHASPAGTYAHTVN RNWYSDADMPASARQEGCQDIATQLISNMDIDVILGGGRKYMFPMPGTPDPEY PADASQNGIRLDGKNLVQEWLAKHQGAWYVWNRTELMQASLDQSVTHLMGL FEPGDTKYEIHRDPTLDPSLMEMTEAALRLLSRNPRGFYLFVEGGRIDHGHHE GVAYQALTEAVMFDDAIERAGQLTSEEDTLTLVTADHSHVFSFGGYTLRGSSIF GLAPSKAQDSKAYTSILYGNPGYVFNNSGVRPDVNESESGSPDYQQAAVPL SSETHGGEDVAVFARGPQAHLVHGVQEESFVAHVMAFAACLEPYTACDLAPP AC
Source	Yeast
Target Names	ALPI
Protein Names	Recommended name: Intestinal-type alkaline phosphatase Short name= IAP Short name= Intestinal alkaline phosphatase EC= 3.1.3.1
Expression Region	26-500aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.