



Recombinant Vaccinia virus Protein E6 (VACWR062), partial

Product Code	CSB-MP325116VAI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21607
Product Type	Recombinant Protein
Immunogen Species	Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strain WR))
Purity	>85% (SDS-PAGE)
Sequence	IWAYLSKKDTGIEFADNDRQDIYTLFQQTGRIVHSNLTETFRDYIFPGDKTSYW VWLNESIANDADIVLNRHAITMYDKILSYIYSEIKQGRVKNMMLKLVYIFEPEKDI RELLLEIIYDIPGDILSIIDAKNDDWKKYFISFYKANFINGNTFISDRTFNEDLFRVV VQIDPEYFDNERIMSLFSTSAADIKRFEDELINNSYISNIIYEVNDITLDTMDDMK KCQIFNEDTSYYVKEYNTYLFLHESD
Source	Mammalian cell
Target Names	VACWR062
Protein Names	Recommended name: Protein E6
Expression Region	180-429aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.