



# Recombinant Mouse Type-2 angiotensin II receptor (Agtr2)

<b>Product Code</b>	CSB-CF001466MO(A4)
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P35374
<b>Alias</b>	Angiotensin II type-2 receptor; AT2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	MKDNFSFAATSRNITSSRPFNDLNATGTNESAFNC SHKPSDKHLEAIPVLYYMI FVIGFA VNIVVVS L FCCQKGP KKVSS IYIFNLALADLLLLATLPLWATYYSYRYDWLFGPV MCKVF GSFLTLNMFASIFFITCMSVDRYQSVIYPFLSQRNPWQASYVVPLVWCMACL SSLPTFY FRDVRTIEYLGVNACIMAFPPEKYAQWSAGIALMKNILGFIPLIFIATCYFGIRKH LLK TNSYGKNRITRDQVLKMAAAVVLAFIICWLPFHVLTFLDALTWMGIINSCEVIAVI DLAL PFAILLGFTNSCVNPFLYCFVGNRFQQKLRSVFRVPITWLQGKRETMSCRKGS SLREMDT FVS
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Agtr2
<b>Protein Names</b>	Type-2 angiotensin II receptor
<b>Expression Region</b>	1-363aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.