




Recombinant Danio rerio D-2-hydroxyglutarate dehydrogenase, mitochondrial (d2hgdh)

Product Code	CSB-EP374600DIL
Relevance	Catalyzes the oxidation of D-2-hydroxyglutarate to alpha-ketoglutarate.
Abbreviation	Recombinant Zebrafish d2hgdh protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1L258
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	ANPANPSAAPPRLPFSRVTQEDLSFFRALLPGRTITDPDLLKSSNVDWLKTVQ GSSDVLLRPKTTEGVSQILRYCNERNLAVCPQGGNTGLVGGSVPVFDEIILSTS LMNQVFAFDNISGILTCQAGCVLENLSHYLEERDFIMPLDLGAKGSCHIGGNVS TNAGGLRLLRYGSLRGTVLGLEVVADGHVNLCLATLRKDNTGYDLKQLFIGS EGTLGVITAVSILCPRKPKAVNVAFLGCSSFQQLLETFQCCRGMLGEILSAFEF LDASCMNLLKHLKLTNPITECPFYVIETAGSNATHDEEKLHQFLEEVMTSSLV TDGTVATEATKIKALWLSLRERVTEALTHEGYTYKYDISLPVEKIYDLVQDMRRH LGGMAKNVVG YGHVGDGNLHLNITSPSKDFDLLAAIEPYVYEWTSQWKGSISA EHGLGLKKNRYIYYSKPSEAVALMGSIKAMLDPKGILNPYKTLDPDNIN
Research Area	Others
Source	E.coli
Target Names	d2hgdh
Expression Region	56-533aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	68.3kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.