



Recombinant Human Rab GDP dissociation inhibitor beta (GDI2), partial

Product Code	CSB-EP009354HU1
Relevance	Regulates the GDP/GTP exchange reaction of most Rab proteins by inhibiting the dissociation of GDP from th, and the subsequent binding of GTP to th.
Abbreviation	Recombinant Human GDI2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50395
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MNEEYDVIVLGTGLTECILSGIMSVNGKKVLHMDRNPYYGGESASITPLEDLK RFKIPGSPPEMGRGRDWNVDLIPKFLMANGQLVKMLLYTEVTRYLDFKVFEG SFVYKGGKIYKVPSTEAELASSLMGLFEKRRFRKFLVYVANFDEKDPRTFEG DPKKTMRDVYKFDLGGQDVIDFTGHALALYRTDDYLDQPCYETINRIKLYSES LARYGKSPYLYPLYGLGELPQGFARLSAIYGGTYMLNKPIEEIIVQNGKVIGVKS EGEIARCKQLICDPSYVKDRVEKVGQVIRVICILSHPIKNTNDANSCQIIIPQNQV NRKSDIYVCMISFAHNVAQAQGYIAIVSTTVETKEPEKEIRPALELLEPIEQKFVSI SDLLVPKDLGTESQIFISRTYDATTHFETTCDDIKNIYKRMTGSEFD FEEMKRKK NDI
Research Area	Signal Transduction
Source	E.coli
Target Names	GDI2
Expression Region	1-441aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	77.2kDa
Protein Length	Partial
Image	

