




Recombinant Human Cathepsin D (CTSD), partial

Product Code	CSB-EP006187HU1
Relevance	Acid protease active in intracellular protein breakdown. Involved in the pathogenesis of several diseases such as breast cancer and possibly Alzheimer disease.
Abbreviation	Recombinant Human CTSD protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07339
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	IPEVLKNYMDAQYYGEIGIGTPPQCFTVVFDTGSSNLWVPSIHCKLLDIACWIH HKYNSDKSSTYVKNGTSFDIHYGSGSLSGYLSQDTVSVPCQSASSASALGGV KVERQVFGKATKQPGITFIAAKFDGILGMAYPRISVNNVLPVFDNLMQQKLVDQ NIFSFYLSRDPDAQPGGELMLGGTDSKYYYKGSLSYLNVTRKAYWQVHLDQVE VASGLTLCKEGCEAIVDTGTSLMVGPDDEVRELQKAIGAVPLIQGEYMIPCEKV STLPAILKLGKGYKLSPEYTLKVSQAGKTLCLSGFMGMDIPPPSGPLWILG DVFIGRYYTVFDRDNNR
Research Area	Neuroscience
Source	E.coli
Target Names	CTSD
Expression Region	67-403aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	40.8kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.