




Recombinant Human COMM domain-containing protein 4 (COMMD4), partial

Product Code	CSB-EP880932HU1
Relevance	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes . Down-regulates activation of NF-kappa-B.
Abbreviation	Recombinant Human COMMD4 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0A8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MRFRFCGDLDCPDWVLAIEISTLAKMSSVKLRLLCSQVLKELLGQGIDYEKILKL TADAKFESGDVKATVAVLSFILSSAAKHSVDGESLSSELQQGLPKEHAASLCR CYEEKQSPLQKHLRVCSLRMNRLAGVGWRVDYTLSSSLLQSVEEPMVHLRLE VAAAPGTPAQPVAMSLSADKFQVLLAELKQAQTLM
Research Area	Cell Biology
Source	E.coli
Target Names	COMMD4
Expression Region	1-195aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	48.4kDa
Protein Length	Partial
Image	 (Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.