




Recombinant Human C-X-C motif chemokine 5 (CXCL5), partial

Product Code	CSB-EP006250HU1
Relevance	Involved in neutrophil activation. In vitro, ENA-78(8-78) and ENA-78(9-78) show a threefold higher chotactic activity for neutrophil granulocytes.
Abbreviation	Recombinant Human CXCL5 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42830
Alias	ENA-78(1-78)Epithelial-derived neutrophil-activating protein 78Neutrophil-activating peptide ENA-78Small-inducible cytokine B5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	AGPAAAVLRELRCVCLQTTQGVHPKMISNLQVFAIGPQCCKVEVVASLKNGKE ICLDPEAPFLKKVIQKILDGG
Source	E.coli
Target Names	CXCL5
Expression Region	37-110aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	11.9kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

