




# Recombinant Mouse C-C motif chemokine 9 (Ccl9), partial

|                          |                                                                                                                                                                                                                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-RP091194m                                                                                                                                                                                                                                                                                    |
| <b>Relevance</b>         | Monokine with inflammatory, pyrogenic and chokinetic properties. Circulates at high concentrations in the blood of healthy animals. Binding to a high-affinity receptor activates calcium release in neutrophils. It also inhibits colony formation of bone marrow myeloid immature progenitors. |
| <b>Abbreviation</b>      | Recombinant Mouse Ccl9 protein, partial                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                  |
| <b>Uniprot No.</b>       | P51670                                                                                                                                                                                                                                                                                           |
| <b>Alias</b>             | CCF18Macrophage inflammatory protein 1-gamma ;MIP-1-gammaMacrophage inflammatory protein-related protein 2 ;MRP-2Small-inducible cytokine A9                                                                                                                                                     |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                             |
| <b>Purity</b>            | ≥ 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                 |
| <b>Sequence</b>          | QITHATETKEVQSSLKAQQGLEIEMFHMGMFQDSSDCCLSYNSRIQCSRFIGYF<br>PTSGGCTRPGIIFISKRGFQVCANPSDRRVQRCIERLEQN                                                                                                                                                                                               |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                                           |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>      | Ccl9                                                                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 22-114aa                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                              |
| <b>Tag Info</b>          | N-terminal 6xHis-tagged                                                                                                                                                                                                                                                                          |
| <b>Mol. Weight</b>       | 14.6kDa                                                                                                                                                                                                                                                                                          |
| <b>Protein Length</b>    | Partial                                                                                                                                                                                                                                                                                          |
| <b>Image</b>             |                                                                                                                                                                                                                                                                                                  |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.