




Recombinant Human Fibroblast growth factor 5 (FGF5),partial

| | |
|--------------------------|--|
| Product Code | CSB-EP008632HU1 |
| Relevance | Plays an important role in the regulation of cell proliferation and cell differentiation. Required for normal regulation of the hair growth cycle. Functions as an inhibitor of hair elongation by promoting progression from anagen, the growth phase of the hair follicle, into catagen the apoptosis-induced regression phase . |
| Abbreviation | Recombinant Human FGF5 protein, partial |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P12034 |
| Product Type | Recombinant Proteins |
| Immunogen Species | Homo sapiens (Human) |
| Purity | Greater than 90% as determined by SDS-PAGE. |
| Sequence | AWAHGEKRLAPKGQPGPAATDRNPRGSSSRQSSSSAMSSSSASSSPAASLG SQGSGLEQSSFQWSPSGRRTGSLYCRVGIGFHLQIYPDGKVNGSHEANMLSV LEIFAVSQGIVGIRGVFSNKFLAMSKKGLHASAKFTDDCKFRERFQENSYNTY ASAIHRTEKTGREWYVALNKRKAKRGCSPRVKPQHISTHFLPRFKQSEQPEL SFTVTVPEKKKPPSPIKPKIPLSAPRKNTNSVKYRLKFRFG |
| Research Area | Cardiovascular |
| Source | E.coli |
| Target Names | FGF5 |
| Expression Region | 18-268aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal GST-tagged |
| Mol. Weight | 54.6kDa |
| Protein Length | Partial |
| Image | |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.