



Recombinant Mouse Renin-2 (Ren2), partial

Product Code CSB-EP360582MO1

Relevance

Adipocyte-secreted protein (adipokine) that regulates adipogenesis, metabolism and inflammation through activation of the chokine-like receptor 1 (CMKLR1). Its other ligands include G protein-coupled receptor 1 (GPR1) and chokine receptor-like 2 (CCRL2). Positively regulates adipocyte differentiation, modulates the expression of adipocyte genes involved in lipid and glucose metabolism and might play a role in angiogenesis, a process essential for the expansion of white adipose tissue. Also acts as a proinflammatory adipokine, causing an increase in secretion of proinflammatory and prodiabetic adipokines, which further impair adipose tissue metabolic function and have negative systemic effects including impaired insulin sensitivity, altered glucose and lipid metabolism, and a decrease in vascular function in other tissues. Can have both pro- and anti-inflammatory properties depending on the modality of enzymatic cleavage by different classes of proteases. Acts as a chotactic factor for leukocyte populations expressing CMKLR1, particularly immature plasmacytoid dendritic cells, but also immature myeloid DCs, macrophages and natural killer cells. Exerts an anti-inflammatory role by preventing TNF/TNFA-induced VCAM1 expression and monocytes adhesion in vascular endothelial cells. The effect is mediated via inhibiting activation of NF-kappa-B and CRK/p38 through stimulation of AKT1/NOS3 signaling and nitric oxide production. Its dual role in inflammation and metabolism might provide a link between chronic inflammation and obesity, as well as obesity-related disorders such as type 2 diabetes and cardiovascular disease. Exhibits an antimicrobial function in the skin

Abbreviation Recombinant Mouse Ren2 protein, partial

Storage The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Uniprot No. P00796

Alias ChemerinRAR-responsive protein TIG2Tazarotene-induced gene 2 protein

Product Type Recombinant Protein

Immunogen Species Mus musculus (Mouse)


Purity ≥ 90% as determined by SDS-PAGE.

Sequence SSLTDLISPVVLTNYLNSQYYGEIGIGTPPQTFKVFVFDTGGSANLWVPSTKCSRLY LACGIHSLYESSDSSSYMENGDDFTIHYGSGRVKGFSLQDSVTVGGITVTQTF GEVTELPLIPFMLAQFDGVLGMGFPAQAVGGVTPVFDHILSQGVLEKVFVSY YNRGPHLLGGEVVLGGSDPEHYQGDFHYVSLSKTDSWQITMKGVSVGSSTLL CEEGCEVVVDTGSSFISAPTSSLKLIMQALGAKEKRLHEYVVSQVPTLPDIS FNLGGRAYTSLSTDYVLQYPN

Research Area Neuroscience



Source	E.coli
Target Names	Ren2
Expression Region	64-351aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	35.0kDa
Protein Length	Partial

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.