




# Recombinant Human Metalloproteinase inhibitor 3 (TIMP3), partial

<b>Product Code</b>	CSB-EP023562HU1
<b>Relevance</b>	Complexes with metalloproteinases (such as collagenases) and irreversibly inactivates th by binding to their catalytic zinc cofactor. May form part of a tissue-specific acute response to remodeling stimuli. Known to act on MMP-1, MMP-2, MMP-3, MMP-7, MMP-9, MMP-13, MMP-14 and MMP-15.
<b>Abbreviation</b>	Recombinant Human TIMP3 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P35625
<b>Alias</b>	Protein MIG-5 Tissue inhibitor of metalloproteinases 3 ; TIMP-3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	HPQDAFCNSDIVIRAKVVGKKLKLVKEGPFGLVYTIKQMKMYRGFTKMPHVQYI HTEASESLCGLKLEVNKYQYLLTGRVYDGMKMYTGLCNFVERWDQLTLSQRKG LNYRYHLGCNCKIKSCYYLPCFVTSKNECLWTDMLSNFGYPGYQSKHYACIR QKGGYCSWYRGWAPPDKSIINA
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	TIMP3
<b>Expression Region</b>	30-208aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	24.8kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.