




Recombinant Escherichia coli Ornithine carbamoyltransferase chain F (argF)

Product Code	CSB-RP129074Ba
Relevance	Reversibly catalyzes the transfer of the carbamoyl group from carbamoyl phosphate (CP) to the N(epsilon) atom of ornithine (ORN) to produce L-citrulline, which is a substrate for argininosuccinate synthetase, the enzyme involved in the final step in arginine biosynthesis.
Abbreviation	Recombinant E.coli argF protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06960
Alias	OTCase-2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SDLYKKHFLKLLDFTPAQFTSLLTLAAQLKADKKNNGKEVQKLTGKNIALIFEKDS TRTRCSFEVAAFDQGARVITYLGPSSQIGHKESIKDTARVLGRMYDGIQYRGH GQEVVETLAQYAGVPVWNGLTNEFHPTQLLADLMTMQEHLPGKAFNEMTLVY AGDARNMGNMSEAAALTGLDLRLLAPKACWPEESLVAECSALAEKHGGKI TLTEDVAAGVKGADFIYTDVWVSMGEAKEKWAERIALLRGYQVNAQMMALTD NPNVKFLHCLPAFHDDQTTLGKQMAKEFDLHGGMEVTDEVFESAASIVFDQA ENRMHTIKAVMMATLGE
Research Area	Others
Source	E.coli
Target Names	argF
Expression Region	2-334aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	40.7kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.