




Recombinant Human Zinc finger and SCAN domain-containing protein 21 (ZSCAN21), partial

Product Code	CSB-RP137874h
Relevance	Strong transcriptional activator.
Abbreviation	Recombinant Human ZSCAN21 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y5A6
Alias	Renal carcinoma antigen NY-REN-21 Zinc finger protein 38 homolog ;Zfp-38
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MMTKVLGMAPVLGPRPPQEQVGPLMVKVEEKEEKGYLPSLEMFRQRFRQF GYHDTPGPREALSQRLVLCCEWLRPEIHTKEQILELLVLEQFLTILPQELQAWV QEHCPESAEEAVTLLEDLERELDEPGHQVSTPPNEQKPVWEKISSSGTAKESP SSMQPQPLETSHKYESWGPLYIQESGEEQEFAQDPRKVRDCRLSTQHEESA DEQKGSEAEGKGDIIISVIIANKPEASLERQCVNLENEKGTKPPLQEAGSKKGR ESVPTKPTPGERRYICAECGKAFSNSSNLTKHRRTHHTGEKPYVCTKCGKAFSH SSNLTLHYRTHLVDRPYDCKCGKAFGQSSDLLKHQRMHTEEAPYQCKDCGK AFSGKGLIRHYRIHTGEKPYQCNECGKSFSQHAGLSSHQRLHTGEKPYKCK ECGKAFNHSSNFNKHHRHTGEKPYWCHHCGKTFCSKSNLSK
Research Area	Transcription
Source	E.coli
Target Names	ZSCAN21
Expression Region	1-461aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	56.4kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.