




# Recombinant Helicobacter pylori Orotidine 5'-phosphate decarboxylase (pyrF)

<b>Product Code</b>	CSB-RP143374Ba
<b>Relevance</b>	Catalyzes the decarboxylation of orotidine 5'-monophosphate (OMP) to uridine 5'-monophosphate (UMP).
<b>Abbreviation</b>	Recombinant Helicobacter pylori pyrF protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P56155
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MQLCVALDLEKKEDNLSLLQELKGLDLWAKVGLRSFIRDGAVFLDEIRKIDENF KIFLDLKLYDIPYTMANAALCAKLDIDMLTVHLSSAKSALTALMQRLNALKKRP LIMGVSALTSFSEEEFLMVYNAPLKTQAIKLSAMGKESGIDGVVCSVFESLAIKE ALGKDFLTLTPGIRLDQNDKEDQERVANAKEAKQNLSDFIVVGRPIYQAKEPRE VVLELLKDC
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	pyrF
<b>Expression Region</b>	1-227aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	29.3 kDa
<b>Protein Length</b>	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.