




Recombinant Salmonella typhi Protein viaA (viaA), partial

Product Code	CSB-RP165194Ba
Abbreviation	Recombinant Salmonella typhi viaA protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8Z2R0
Alias	VWA domain protein interacting with AAA ATPase
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhi
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	TLDTLNTMLAVSEEGMVEEMILALLASPQLVIFFEKFPRLKNAVTADLPRWREA LRSRLKDARVPPELTEEVMCYQQSQLLSTPQFIVQLPQILALLHRLHSPYAAQA KQLTESNSTFTPALHTLFLQRWRLSLVVQATTLNQQLLEEEREQLLSDVQERM TLSGQLEPTLAENDNAAGRLWDMSAGQLKRGDYQLIVKYGEFLAAQPELMQL AEQLGRSREAKSVPKKDAPMETFRTLVPATVPEQVDGIQQGDDILRLPPE LATLGITELEYEFYRRLVEKQLLTYRLHGEAWREKVTERPVVHQDVDEQPRGP FIVCVDTS GSMGGFNEQCAKAFCLALMRVALADNRRCFIMLFSTDVVRYELSG PESIEQAIRFLSQRFRRGGTDIASCFAIIRMQGREWFDADAVVISDFIAQRLPD DVVSKVGELQRLHQHFRHAVAMSAHGKPGIMRIFDHIWRFDTGMRSR
Research Area	Others
Source	E.coli
Target Names	viaA
Expression Region	3-476aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	58.2kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.