




# Recombinant *Thermus aquaticus* DNA polymerase I, thermostable (polA), partial

<b>Product Code</b>	CSB-EP323156TNN4
<b>Abbreviation</b>	Recombinant <i>Thermus aquaticus</i> polA protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P19821
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	<i>Thermus aquaticus</i>
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	LLESPKALEEAPWPPPEGAFVGFVLSRKEPMWADLLALAAARGGRVHRAPEP YKALRDLKEARGLLAKDLSVLALREGLGLPPGDDPMLLAYLLDPSNTTPEGVA RRYGGEWTEEAGERAALSERLFANLWGRLEGEERLLWLYREVERPLSAVLAH MEATGVRLDVAYLRALSLEVAEEIARLEAEVFRLAGHPFNLNSRDQLERVLDFE LGLPAIGKTEKTGKRSTSAAVLEALREAHPIVEKILQYRELTCLKSTYIDPLPLDI HPRTGRLHTRFNQTATATGRLSSSDPNLQNPVVRTPLGQRIRRAFIAEEGWLLV ALDYSQIELRVLAHLSGDENLIRVFQEGRDIHTETASWMFGVPREAVDPLMRR AAKTINFGVLYGMSAHRLSQELAIPIYEEAQAFIERYFQSFPKVRWIEKTLLEEG RRRGYVETLFGRRRYVPDLEARVKSVDREAAERMAFNMPVQGTAAADLMKLAM VKLFPRLEEMGARMMLLQVHDELVLEAPKERAEAVARLAKEVMEGVYPLAVPLE VEVGIGEDWLSAKE
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	polA
<b>Expression Region</b>	287-832aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	88.5kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.