




Recombinant *Saccharomyces cerevisiae* Neutral trehalase (NTH1), partial

Product Code	CSB-EP330157SVG1
Abbreviation	Recombinant <i>Saccharomyces cerevisiae</i> NTH1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32356
Product Type	Recombinant Proteins
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QVNTSQGPVAQGRQRRLLSSLSEFNDFPFSNAEVYYPPTDPRKQKQAKPAKIN RTRTMSVFDNVSPFKKTGFGKLQQTRRGSEDDTYSSSQGNRRFFIEDVDKTL NELLAEDTDKNYQITIEDTGPKVLKVGTTANSYGYKHINIRGTYMLSNLLQELTI AKSFGRHQIFLDEARINENPVNRLSRLINTQFWSLTRRVDLNNVGEIAKDTKID TPGAK
Research Area	Others
Source	E.coli
Target Names	NTH1
Expression Region	3-221aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	28.8kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.