




Recombinant Hepatitis C Genome polyprotein

Product Code	CSB-RP184674Ba
Abbreviation	Recombinant Hepatitis C Genome polyprotein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	S4UAW6
Product Type	Recombinant Protein
Immunogen Species	Hepatitis C
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	APITAYAQQTRGLFSTIVTSLTGRDTNENRGEVQVLSTATQSFLGTAVNGVMW TVYHGAGSKTISGPKGPVNQMYTNVDQDLVGWPAPPGVRSLAPCTCGASDL YLVTRHADVVPVRRRGDTRGALLSPRIATLKGSSGGPLLCPMGHAAGLFRAA VCTRGVAKAVDFVSVESLETTMRSPVFTDNSTPPAVPQTYQVAHLHAPTGSG KSTKVPAAAYAAQGYKVLVLNPSVAATLGFGAYMSKAYGIDPNIRTGVRTITTGA PITYSTYGKFLADGGCSGGAYDIIICDECHSTDSTTILGIGTVLDQAETAGVRLTI LATATPPGSVTTPHSNIEEVALPTTGEIPFYGKAIPLELIKGGRLIFCHSKKKCD ELAKQLTSLGLNAVAYYRGLDVSVIPTSGDVVVCATDALMTGYTGDFDSVIDC NTSVIQTVDFSLDPTFSIETTTVPQDAVSRSQRRGRTGRGRLGIYRYVTPGER PSGIFDTSVLCECYDAGCAWYELTPAETTTRLRAYFNTPGLPVCQDHLEFWES VFTGLTQIDGHFLSQTKQTGENFPYLVAYQATVCARALAPPPSWDTMWKCLIR LKPTLHGPTPLLYRLGSVQNEVTLTHPITKYIMACMSADLEVVT
Research Area	Others
Source	E.coli
Expression Region	1-631aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	71.3 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.