




Recombinant Escherichia coli 4-hydroxy-tetrahydrodipicolinate synthase (dapA)

Product Code	CSB-EP363880ENV
Relevance	Catalyzes the condensation of (S)-aspartate-beta-sialdehyde [(S)-ASA] and pyruvate to 4-hydroxy-tetrahydrodipicolinate (HTPA).
Abbreviation	Recombinant E.coli dapA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A6L2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MFTGSIVAIVTPMDEKGNVCRASLKKLIDYHVASGTSIAIVSVGTTGESATLNHD EHADVMMMLDLADGRIPVIAGTGANATAEAIISLTQRFNDSGIVGCLTVTPYYN RPSQEGLYQHFKAIAEHTDLPQILYNVPSRTGCDLLPETVGRRLAKVKNIIGIKEA TGNLTRVNQIKELVSDDFVLLSGDDASALDFMQLGGHGVISVTANVAARDMAQ MCKLAAEGHFAEARVINQRLMPLHNKLFVEPNPIPVKWACKELGLVATDTLRL PMTPIITDSGRETVRAALKHAGLL
Research Area	Others
Source	E.coli
Target Names	dapA
Expression Region	1-292aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	35.3kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.