




Recombinant Human Paraneoplastic antigen Ma2 (PNMA2), partial

Product Code	CSB-YP018267HU
Abbreviation	Recombinant Human PNMA2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UL42
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ALALLEDWCRIMSVDEQKSLMVTGIPADFEEAEIQEVLQETLKSLGRYRLLGKIF RKQENANAVLLELLEDTDVSAIPSEVQKGKGVWKVIFKTPNQDTEFLERLNLFL EKEGQTVSGMFRALGQEGVSPATVPCISPELLAHLGQAMAHAPQPLPMRY RKLRVFGSAVPAPEEESFEVWLEQATEIVKEWPVTEAEKKRWLAESLRGPAL DLMHIVQADNPSISVEECLEAFKQVFGSLESRRTAQVRYLKTYQEEGEKVSAY VLRLETLLRRAVEKRAIPRRIADQVRLEQVMAGATLNQMLWCRLRELKDQGPP PSFLELMKVIREEEEEEEASFENESIEEPEERDGYGRWNHE
Research Area	Neuroscience
Source	Yeast
Target Names	PNMA2
Protein Names	Recommended name: Paraneoplastic antigen Ma2Alternative name(s): 40 kDa neuronal protein Onconeural antigen Ma2 Paraneoplastic neuronal antigen MM2
Expression Region	2-261aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	43.1kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.