




# Recombinant Helicobacter pylori Urease subunit alpha (ureA)

<b>Product Code</b>	CSB-YP320452HUV
<b>Relevance</b>	Ammonia produced by ureolysis increases the gastric pH thereby providing an environment permissive for colonization of the stomach.
<b>Abbreviation</b>	Recombinant Helicobacter pylori ureA protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P14916
<b>Alias</b>	Urea amidohydrolase subunit alpha
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MKLTPKELDKLMLHYAGELAKKRKEKGIKLNYVEAVALISAHIMEEARAGKKTA AELMQEGRTLLKPDDVMDGVASMIHEVGVIEAMFPDGTKLVTVHTPIEANGKLV PGELFLKNEDITINEGKKAVSVKVNVDGDRPVQIGSHFFFEVNRCLDFDREKT FGKRLDIASGTAVRFEPGEEKSVELIDIGGNRRIFGFNALVDRQADNESKKIALH RAKERGFHGAKSDDNYVKTIKE
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	ureA
<b>Protein Names</b>	Recommended name: Urease subunit alpha EC= 3.5.1.5 Alternative name(s): Urea amidohydrolase subunit alpha
<b>Expression Region</b>	1-238aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	28.5kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.