




# Recombinant Epstein-Barr virus Envelope glycoprotein B (gB), partial

<b>Product Code</b>	CSB-EP355975EFA
<b>Relevance</b>	Envelope glycoprotein that forms spikes at the surface of virion envelope. Involved in fusion of viral and cellular membranes leading to virus entry into the host cell. Following initial binding of gp350/220 to its receptor, membrane fusion is mediated by the fusion machinery composed at least of gB and the heterodimer gH/gL. May also be involved in the fusion between the virion envelope and the outer nuclear membrane during virion morphogenesis
<b>Abbreviation</b>	Recombinant Epstein-Barr virus gB protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P03188
<b>Alias</b>	GP115 Glycoprotein GP110
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	QTPEQPAPPATTVQPTATRQQTSFPFRVCELSSHGDLFRFSSDIQCPSFGTRE NHTEGLLMVFKDNIIPYSFKVRSYTKIVTNILYNGWYADSVTNRHEEKFSVDSY ETDQMDTIYQCYNVAVKMTKDGLTRVYVDRDGVNITVNLKPTGGLANGVRRYA SQTELYDAPGWLIIWYRTRTTVNCLITDMMAKSNSPFDFVTTTGGQTVEMSPF YDGKKNKETFHERADSFHVRTNYKIV
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	gB
<b>Expression Region</b>	23-260aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	43.3kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.