




Recombinant Cricetulus griseus C-C motif chemokine (LOC100763833)

Product Code	CSB-RP2059DXU
Abbreviation	Recombinant Cricetulus griseus NEU2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	G3GTT2
Product Type	Recombinant Proteins
Immunogen Species	Cricetulus griseus
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QPDVNSPLTCCYSFTA KRIPEKRLESYKRITSSKCPKEAVIFITKLKREICADPK QDWVQTYTKKLDQSQAKSEAATVYKTAPLNANLTHESAVNASTTAFPTTDLRT SVRVTSMTVN
Research Area	Others
Source	E.coli
Target Names	LOC100763833
Expression Region	25-143aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	29.3kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.