



Recombinant Human Protein cereblon (CRBN)


Product Code	CSB-CF842761HU
Relevance	Substrate recognition component of a DCX (DDB1-CUL4-X-box) E3 protein ligase complex that mediates the ubiquitination and subsequent proteasomal degradation of target proteins, such as MEIS2. Normal degradation of key regulatory proteins is required for normal limb outgrowth and expression of the fibroblast growth factor FGF8. May play a role in memory and learning by regulating the assembly and neuronal surface expression of large-conductance calcium-activated potassium channels in brain regions involved in memory and learning via its interaction with KCNT1. Binding of pomalidomide and other thalidomide-related drugs changes the substrate specificity of the human protein, leading to decreased degradation of MEIS2 and other target proteins and increased degradation of MYC, IRF4, IKZF1 and IKZF3
Abbreviation	Recombinant Human CRBN protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96SW2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens(Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MAGEGDQQDAAHNMGNHLPLLPAAESEEDEMEVEDQDSKEAKKPNIINFDTSLPTSHTYLGADMEEFHGRTLHDDDDSCQVIPVLPQVMMILIPGQTLPLQLFHPQEVSMVRNLIQKDRTFAYLAYSNVQEREAQFGTTAEIYAYREEQDFGIEIVKVKAIQRQRFKVLRLRTQSDGIQQAQVQILPECVLPSTMSAVQLESLNKCQIFPSKPVSR EDQCSYKWWQKYQKRKFHCANLT SWPRWLYSLYDAETLMDRIKKQLREWDE NLKDDSLPSNPIDFSYRVAACLPIDDLRIQLLKIGSAIQLRLCELDIMNKCTSLC CKQCQETEITTKNEIFSLSLCGPMAAYVNPHGYVHETLTVYKACNLNLIGRPST EHSWFPGYAWTVAQCKICASHIGWKFTATKKDMSPQKFWGLTRSALLPTIPDTEDEISPDKVILCL
Research Area	Neuroscience
Source	in vitro E.coli expression system
Target Names	CRBN
Expression Region	1-442aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	54.5kDa




Protein Length

Full Length


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-CF842761HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) CRBN.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-CF842761HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) CRBN.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.