




Recombinant Dog Cationic trypsin

Product Code	CSB-CF361971DO
Abbreviation	Recombinant Dog Cationic trypsin protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06871
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	IVGGYTCSRNSVPYQVSLNSGYHFCGGSLINSQWVVSAAHCYKSRIQVRLGE YNIAVSEGGEQFINAAKIIRHPRYNANTIDNDIMLIKLSPPATLNSRVSAIALPKSC PAAGTQCLISGWGNTQSIGQNYPDVLQCLKAPILSDSVCRNAYPGQISSNMMC LGYMEGGKDSCQGDSGGPVVCNGELQGVVSWGAGCAQKKGKPGVSPKVCK YVSWIQQTIAAN
Research Area	others
Source	in vitro E.coli expression system
Expression Region	24-246aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	39.6kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.