



Recombinant Human 5-hydroxytryptamine receptor 1F (HTR1F), partial

| | |
|--------------------------|---|
| Product Code | CSB-MP010886HU1 |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P30939 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | YKIYRAAKTLYHKRQASRIAKEEVNGQVLLESGEKSTKSVSTSYVLEKSLSDPS TDFDKI HSTVRSRLRSEFKHEKSWRRQKISGTREK |
| Source | Mammalian cell |
| Target Names | HTR1F |
| Protein Names | Recommended name: 5-hydroxytryptamine receptor 1F Short name= 5-HT-1F Short name= 5-HT1F Alternative name(s): Serotonin receptor 1F |
| Expression Region | 202-290aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |