



ACYP2 Antibody, HRP conjugated

Product Code	CSB-PA001266LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14621
Immunogen	Recombinant Human Acylphosphatase-2 protein (28-87AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acylphosphatase-2 (EC 3.6.1.7) (Acylphosphatase, muscle type isozyme) (Acylphosphate phosphohydrolase 2), ACYP2, ACYP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ACYP2