



ADRA2A Antibody, FITC conjugated

Product Code	CSB-PA001388LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08913
Immunogen	Recombinant Human Alpha-2A adrenergic receptor protein (218-374AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2A adrenergic receptor (Alpha-2 adrenergic receptor subtype C10) (Alpha-2A adrenoreceptor) (Alpha-2A adrenoceptor) (Alpha-2AAR), ADRA2A, ADRA2R ADRAR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ADRA2A