



ANXA4 Antibody, HRP conjugated

Product Code	CSB-PA001845LB01DO
Abbreviation	Annexin A4
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50994
Immunogen	Recombinant Dog Annexin A4 protein (2-319AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Relevance	Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A4 (36 kDa zymogen granule membrane-associated protein) (ZAP36) (Annexin IV) (Annexin-4) (Lipocortin IV), ANXA4, ANX4
Species	Dog
Research Area	Others
Target Names	ANXA4