



BPI Antibody, FITC conjugated

Product Code	CSB-PA002783LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17213
Immunogen	Recombinant Human Bactericidal permeability-increasing protein (149-362AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bactericidal permeability-increasing protein (BPI) (CAP 57), BPI
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	BPI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.