



CBR1 Antibody, FITC conjugated

Product Code	CSB-PA004586LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16152
Immunogen	Recombinant Human Carbonyl reductase [NADPH] 1 protein (176-268AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonyl reductase [NADPH] 1 (EC 1.1.1.184) (15-hydroxyprostaglandin dehydrogenase [NADP(+)] (EC 1.1.1.197) (NADPH-dependent carbonyl reductase 1) (Prostaglandin 9-ketoreductase) (Prostaglandin-E(2) 9-reductase) (EC 1.1.1.189) (Short chain dehydrogenase/reductase family 21C member 1), CBR1, CBR CRN SDR21C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CBR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.