



CERK Antibody, FITC conjugated

Product Code	CSB-PA005256LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCT0
Immunogen	Recombinant Human Ceramide kinase protein (212-478AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ceramide kinase (hCERK) (EC 2.7.1.138) (Acylsphingosine kinase) (Lipid kinase 4) (LK4), CERK, KIAA1646
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CERK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.