



CYP3A4 Antibody, Biotin conjugated

Product Code	CSB-PA006439LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08684
Immunogen	Recombinant Human Cytochrome P450 3A4 protein (259-441AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 3A4 (EC 1.14.14.-) (1,8-cineole 2-exo-monooxygenase) (EC 1.14.14.56) (Albendazole monooxygenase) (EC 1.14.14.-) (Albendazole sulfoxidase) (CYP11A3) (CYP11A4) (Cholesterol 25-hydroxylase) (EC 1.14.14.-) (Cytochrome P450 3A3) (Cytochrome P450 HLp) (Cytochrome P450 NF-25) (Cytochrome P450-PCN1) (Nifedipine oxidase) (Quinine 3-monooxygenase) (EC 1.14.14.55) (Taurochenodeoxycholate 6-alpha-hydroxylase) (EC 1.14.14.-), CYP3A4, CYP3A3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP3A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.