





DAZL Antibody

Product Code	CSB-PA006514LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92904
Immunogen	Recombinant Human Deleted in azoospermia-like protein (130-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deleted in azoospermia-like (DAZ homolog) (DAZ-like autosomal) (Deleted in azoospermia-like 1) (SPGY-like-autosomal), DAZL, DAZH DAZL1 DAZLA SPGYLA
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DAZL

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA006514LA01HU at dilution of 1:100



Immunofluorescent analysis of U251 cells using CSB-PA006514LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.