



DAZL Antibody, FITC conjugated

Product Code	CSB-PA006514LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92904
Immunogen	Recombinant Human Deleted in azoospermia-like protein (130-250AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deleted in azoospermia-like (DAZ homolog) (DAZ-like autosomal) (Deleted in azoospermia-like 1) (SPGY-like-autosomal), DAZL, DAZH DAZL1 DAZLA SPGYLA
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DAZL