



EPHA4 Antibody, FITC conjugated

Product Code	CSB-PA007724LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54764
Immunogen	Recombinant Human Ephrin type-A receptor 4 protein (345-545AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin type-A receptor 4 (EC 2.7.10.1) (EPH-like kinase 8) (EK8) (hEK8) (Tyrosine-protein kinase TYRO1) (Tyrosine-protein kinase receptor SEK), EPHA4, HEK8 SEK TYRO1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EPHA4