



ERAL1 Antibody, Biotin conjugated

Product Code	CSB-PA007759LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75616
Immunogen	Recombinant Human GTPase Era, mitochondrial protein (271-382AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase Era, mitochondrial (H-ERA) (hERA) (Conserved ERA-like GTPase) (CEGA) (ERA-W) (ERA-like protein 1), ERAL1, HERA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERAL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.