




FAM3A Antibody

Product Code	CSB-PA008226LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98173
Immunogen	Recombinant Human Protein FAM3A protein (151-216AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM3A (Cytokine-like protein 2-19), FAM3A, 2-19 2.19
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	FAM3A

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA008226LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.