



FAM3A Antibody, FITC conjugated

Product Code	CSB-PA008226LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98173
Immunogen	Recombinant Human Protein FAM3A protein (151-216AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM3A (Cytokine-like protein 2-19), FAM3A, 2-19 2.19
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	FAM3A