




Gba Antibody

Product Code	CSB-PA009289HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17439
Immunogen	Recombinant Mouse Glucosylceramidase protein (20-515AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucosylceramidase (EC 3.2.1.45) (Acid beta-glucosidase) (Beta-glucocerebrosidase) (D-glucosyl-N-acylsphingosine glucohydrolase), Gba
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Gba

Image



Western Blot

Positive WB detected in: Mouse liver tissue, Mouse spleen tissue, Mouse lung tissue

All lanes: Gba antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 58 kDa

Observed band size: 58 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.