





GBP2 Antibody

Product Code	CSB-PA009298LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32456
Immunogen	Recombinant Human Guanylate-binding protein 2 protein (100-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC, IF; Recommended dilution:IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	GBP2

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA009298LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA009298LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)