



GBP2 Antibody, Biotin conjugated

Product Code	CSB-PA009298LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32456
Immunogen	Recombinant Human Guanylate-binding protein 2 protein (100-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanylate-binding protein 2 (EC 3.6.5.-) (GTP-binding protein 2) (GBP-2) (HuGBP-2) (Guanine nucleotide-binding protein 2) (Interferon-induced guanylate-binding protein 2), GBP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	GBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.