



GFRA1 Antibody, HRP conjugated

Product Code	CSB-PA009379LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56159
Immunogen	Recombinant Human GDNF family receptor alpha-1 protein (230-385AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDNF family receptor alpha-1 (GDNF receptor alpha-1) (GDNFR-alpha-1) (GFR-alpha-1) (RET ligand 1) (TGF-beta-related neurotrophic factor receptor 1), GFRA1, GDNFRA RETL1 TRNR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GFRA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.