



GFRA1 Antibody, Biotin conjugated

Product Code	CSB-PA009379LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56159
Immunogen	Recombinant Human GDNF family receptor alpha-1 protein (230-385AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDNF family receptor alpha-1 (GDNF receptor alpha-1) (GDNFR-alpha-1) (GFR-alpha-1) (RET ligand 1) (TGF-beta-related neurotrophic factor receptor 1), GFRA1, GDNFRA RETL1 TRNR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GFRA1