



GFRA2 Antibody, Biotin conjugated

Product Code	CSB-PA009380LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00451
Immunogen	Recombinant Human GDNF family receptor alpha-2 protein (242-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDNF family receptor alpha-2 (GDNF receptor alpha-2) (GDNFR-alpha-2) (GFR-alpha-2) (GDNF receptor beta) (GDNFR-beta) (Neurturin receptor alpha) (NRTNR-alpha) (NTNR-alpha) (RET ligand 2) (TGF-beta-related neurotrophic factor receptor 2), GFRA2, GDNFRB RETL2 TRNR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GFRA2